

UV/VIS Spectrophotometer Calibration Standards

Chromate calibration solution		250 ml	CL02.0328.0250
100 mg K ₂ Cr ₂ O ₇ / l H ₂ SO ₄ 0.005 mol/l			
Chromate calibration solution (UV-VIS Standard 1A)	NEW	1 ml	CL02.0377.0001
600.6 mg K ₂ Cr ₂ O ₇ / l H ₂ SO ₄ 0.01 N			
Chromate calibration solution (UV-VIS Standard 1)	NEW	10 ml	CL02.0374.0010
60 mg K ₂ Cr ₂ O ₇ / l H ₂ SO ₄ 0.01 N			
Potassium chloride solution (UV-VIS Standard 4)	NEW	10 ml	CL02.1198.0010
12 g KCl / l in H ₂ O			
Toluene standard solution (UV-VIS Standard 5)	NEW	10 ml	CL41.2065.0010
Solution contains 0,02% (v/v) in n-Hexane			
Holmium(III) oxide calibration solution (UV-VIS Standard 6)	NEW	100 ml	CL02.0804.0100
5% Ho ₂ O ₃ / 1.4 mol/l HClO ₄			