

1

Electrophoresis Buffers, Ready-to-Use

SERVA's ready-to-use electrophoresis buffers save not only time but guarantee best results because they are made from high-quality reagents and application tested.

Product	Size	Cat. no.
TAE Buffer, 10x	1 L	42553.01
	10 L	42553.04
TAE Buffer, 50x, molecular biology grade	1 L	42549.01
TBE Buffer, 10x	2x 500 ml	42557.01

2

Reagents

Reagents "electrophoresis grade" are in-house tested for electrophoresis applications. "Molecular biology grade" reagents are guaranteed DNase/RNase-free.

Product	Size	Cat. no.
Acetic acid 100 %, analytical grade	1 L	45633.01
Boric acid, electrophoresis grade	1 kg	15166.02
	5 g	15375.01
	25 g	15375.02
Bromophenol Blue-Na-salt	100 g	11280.01
	1 kg	11280.02
	5 kg	11280.03
Ethylenediamine tetraacetic acid-Na ₂ -salt, analytical grade	250 g	39760.01
Glycerol from plant 87 % molecular biology grade	1 L	39788.01
	500 g	37186.02
Tris(hydroxymethyl)aminomethane, molecular biology grade	1 kg	37186.03
	2.5 kg	37186.04
	500 g	37181.01
Tris(hydroxymethyl)aminomethane, electrophoresis grade	1 kg	37181.02
	2.5 kg	37181.03
Sucrose, analytical grade	500 g	35579.02
	5 kg	35579.03

Ready-to-use solutions for saving time and work

Application tested for best results in electrophoresis